

Hong ze

POWDER CONTENT TESTING RECORD

Document No.	BA-R-015	Revision Date	2013.08.01	Revision Times	2	Issued Date	2002.12.1
--------------	----------	---------------	------------	----------------	---	-------------	-----------

Product Name: Nitrile Pre-powdered Powder Free Color: blue
 Textured Smooth Others: /

Customer Name/Production Line: SETINO PO No./ Lot Number: 01-G-3

Testing Date: 2020.5.24 Testing Time: 10:00

Testing Equipment #: Incubator: HZ3C-013-002 Electronic balance: HZ3C-015-001

Test Record: Take the average result of five gloves,
 Unit: For pre-powdered glove, the result for powder content should be converted to mg/dm².
 For Powder free glove, the unit for powder content is mg/glove.

Size:	<u>M</u>	<u>L</u>				
Average content per glove	<u>0.54mg</u>	<u>0.84mg</u>				
Average content per dm ²	<u>/</u>	<u>/</u>				

Powder Content Criteria: Not more than 10mg/dm² for pre-powdered glove.

Not more than 2mg/glove for powder free glove.

Test Result:



Fail

QA Director: Xu Li Hua

Inspector: DI NA