

# TEST REPORT: 7191240304-CHM20/01-JS

Date: 17 JUL 2020

Tel: +65 68851312 Fax: +65 67784301

Client's Ref: 221234218

Email: zhou.xiao@tuv-sud-psb.sg

**Note:** This report is issued subject to the Testing and Certification Regulations of the TÜV SÜD Group and the General Terms and Conditions of Business of TÜV SÜD PSB Pte Ltd. In addition, this report is governed by the terms set out within this report.



PSB Singapore

Add value.  
Inspire trust.

## SUBJECT

Determination of Removable Surface Powder for Glove Samples.

## CLIENT

Sri Trang Gloves (Thailand) Public Company Limited  
110 Kanchanawanit., T.Patong, A.Hatyai,  
Songkhla 90230  
Thailand

Attention: Miss Nattatida Promsila

## SAMPLE SUBMISSION DATE

03 Jul 2020

## DESCRIPTION OF SAMPLE

One glove sample labelled as follows was received.

S/N	Product Description	Colour	Size
1	Latex Powder Free Online Chlorination Gloves (013145956)	Natural	M

## DATE OF ANALYSIS

08 Jul 2020 – 16 Jul 2020

## METHOD OF TEST

The powder residual test was performed in accordance *EN 455-3:2015 Medical glove for single use Part 3: Requirements and testing for biological evaluation, clause 5.2 powder* by using test method *ISO 21171:2006: Medical gloves – Determination of removable surface powder, Clause 9.2. Method B - Procedure for "powder-free" gloves other than surgeon's gloves.*



TÜV SÜD PSB

Laboratory:  
TÜV SÜD PSB Pte. Ltd.  
No.1 Science Park Drive  
Singapore 118221

Phone : +65-6885 1333  
Fax : +65-6776 8670  
E-mail: enquiries@tuv-sud-psb.sg  
www.tuv-sud-psb.sg  
Co. Reg : 199002667R

Regional Head Office:  
TÜV SÜD Asia Pacific Pte. Ltd.  
1 Science Park Drive, #02-01  
Singapore 118221  
**TÜV®**

# TEST REPORT: 7191240304-CHM20/01-JS

17 JUL 2020



PSB Singapore

## RESULTS

**Table 1: Analytical Result of Removable Surface Powder for Latex Powder Free Online Chlorination Gloves Sample.**

Sample		Test Result	EN 455-3 Requirement	Inferred Result
1	Latex Powder Free Online Chlorination Gloves (013145956)		Powder residue limit $\leq 2.0$ mg/glove for Powder-Free Gloves	Pass
	Colour	Natural		
	Size	M		

**MS SIM SIEW HEONG**  
TECHNICAL EXECUTIVE

**DR XIAO ZHOU**  
PRODUCT MANAGER  
MICROCONTAMINATION DIAGNOSIS  
CHEMICAL & MATERIALS



# TEST REPORT: 7191240304-CHM20/01-JS

17 JUL 2020



PSB Singapore

Please note that this Report is issued under the following terms :

1. This report applies to the sample of the specific product/equipment given at the time of its testing/calibration. The results are not used to indicate or imply that they are applicable to other similar items. In addition, such results must not be used to indicate or imply that TÜV SÜD PSB approves, recommends or endorses the manufacturer, supplier or user of such product/equipment, or that TÜV SÜD PSB in any way "guarantees" the later performance of the product/equipment. Unless otherwise stated in this report, no tests were conducted to determine long term effects of using the specific product/equipment.
2. The sample/s mentioned in this report is/are submitted/supplied/manufactured by the Client. TÜV SÜD PSB therefore assumes no responsibility for the accuracy of information on the brand name, model number, origin of manufacture, consignment or any information supplied.
3. Nothing in this report shall be interpreted to mean that TÜV SÜD PSB has verified or ascertained any endorsement or marks from any other testing authority or bodies that may be found on that sample.
4. This report shall not be reproduced wholly or in parts and no reference shall be made by the Client to TÜV SÜD PSB or to the report or results furnished by TÜV SÜD PSB in any advertisements or sales promotion.
5. Unless otherwise stated, the tests were carried out in TÜV SÜD PSB Pte Ltd, No.1 Science Park Drive Singapore 118221.

July 2011

